

# Protocols

Krishna  
17-11-99

- Slide treatment with pepsin
- DNA extraction from YACs
- PAC and BAC minipreps
- ALU-PCR amplification
- Probe labelling
- Fluorescence *in situ* hybridization
- References

## Slide treatment with pepsin (for ageing fresh made slides)▲

- incubate for 30' at 37°C in '0.005% pepsin/0.01M HCl'\* (stock)
- wash 1xPBS at RT for 5'
- 5' at RT in buffer\* for 50ml\*: 5ml 10xPBS  
5ml 0.5M MgCl<sub>2</sub>  
40ml H<sub>2</sub>O)
- 5' in paraformaldehyde 4% at RT\* for 50ml: 5ml 10xPBS  
5ml 0.5M MgCl<sub>2</sub>  
15ml H<sub>2</sub>O  
25ml paraformaldehyde 8%
- wash 1xPBS, 5' at RT
- 3 wash in ethanol series (70%, 90%, and 100%)
- air drying

### Rapid ageing of slides:

- prepare slides as usual
- ON incubation at 37°C
- 2h in oven at 80-90°C
- pepsin treatment

---

## DNA extraction from YACs (CEPH megalibrary)▲

- seed the YAC in 5-10ml of YPD (with ampicillin)
- incubate 2 days at 30°C in a stirrer
- fuge 5' at 3500RPM; discard the SN
- suspend in 0.5ml 1M sorbitol-0.1M EDTA (pH7.5);
- transfer in eppendorf tube
- add 7ul of lyticase (stock solution: 10.000U/ml).
- incubate at 37°C for 30'
- fuge for 1', discard the SN.
- suspend in 0.5ml of 50mM Tris-HCl (pH7.4) - 20mM EDTA
- add 50ul of SDS 10%; mix

- at 65°C for 30'
- add 0.2ml of 5M KAc; in ice for 60'
- fuge for 5'
- transfer the SN in a eppendorf tube
- add 1 vol isopropanol at RT
- mix; leave at RT per 5'
- fuge for 10 seconds, discard the SN, let the pellet dry
- suspend in 0.3ml of TE (pH7.4)
- add 15µl of RNase A 1mg/ml; at 37°C for 30'
- add 30µl of Na Acet., mix
- precipitate with 0.2ml of isopropanol
- fuge briefly; discard the supernatant
- dry the pellet in Savant for 5'
- suspend in 0.1-0.3ml TE (pH7.4)
- fuge briefly to eliminate insoluble material

-----  
**YPD:** (1 liter):     10gr yeast extract  
                   &nbsp; nbsp;       &nbsp; nbsp;     20gr peptone  
 -----

Autoclave for 30'; when the temperature is around 50C add  
 50ml of glucose (40%, sterilized by filtration)

---

## PAC-BAC Miniprep ▲

**NOTE:** Resistance: BAC: Cloramphenicol; PAC: Kanamycin

1. Grow bacteria in 10 ml bacterial medium containing the right antibiotic in 50 ml falcon tubes for 16/20 hours.
2. Centrifuge 10' at 4000 rpm.
3. Add 300µl of GTE (Glucose 50mM, Tris pH=8 25mM, EDTA 10Mm)
4. Resuspend pellet completely.
5. Trasfer the cell suspension into 2.2ml eppendorf tubes.
6. Add 600 µl of denaturation solution freshly made; do not vortex, mix by inverting several times, do not lyse for more than 5'. The lysate should appear viscous.
7. Add 500µl of 7.5M Ammonium Acetate; do not vortex, mix immediately by inverting several times; leave on ice for 10'; invert several times during the incubation period.
8. Centrifuge at 13000 rpm for 20'
9. Pour the supernatant into fresh 2.2ml eppendorf tubes; the supernatant is most often not clear and a second centrifuge at 14000 rpm for 10' needs to be carried out.
10. Pour the supernatant into fresh 2.2ml eppendorf tubes.
11. Add 700µl isopropanol; mix by inverting several times.
12. Centrifuge at 14000 rpm for 20'.
13. Discard supernatant; the pellet should be barely visible
14. Wash the pellet with 500µl of 70% EtOH.
15. Centrifuge at 14000 rpm for 5'
16. Discard the supernatant; don't let the pellet dry
17. Resuspend the pellet in 100µl of TE by tapping the tubes
18. Treat with RNase f.c. 100ug/ml: at 37C for 30'

19. Precipitate the DNA with 1/10 vol NaAc and 3 vol Ethanol
20. Incubate at -20C for 20'
21. Fuge for 15' at 14000RPM; wash with 70% Ethanol
22. Resuspend in TE (appropriate volume)
23. Check on gel; store at 4 degrees

-----  
Before labeling: fuge for 15' at 14.000RPM

---

## Alu-PCR amplification

PCR amplifications are performed using the following Alu primers (Liu et al., 1993):

5' GGATT ACAGG YRTGA GCCA 3';  
5' RCCAY TGCAC TCCAG CCTG 3' (Y=C/T; R=A/G),

50-100ng of genomic DNA from the hybrid is amplified in a volume of 50ul. Cycles: 5 min at 95°C; then 30 cycles as follows: 1 min at 95°C, 1 min at 65°C, 4 min at 72°C, with a final extension of 10 min at 72°C.

Check the concentration of amplified products on gel. A good amplification yields about 200ng/ul

Label the products using Nick-translation as described (see below).

Reagents for PCR (store frozen):

- 10x dNTPs mix 2mM each, pH7.
- 10x reaction buffer (usually comes with the Taq).

---

## Probe labelling with Biotin by nick-translation

1. Add to a microfuge tube, on ice:
  - 2µg of amplification products (not purified)
  - 10µl 10x nick translation buffer\*
  - 10µl dNTPs mix\*
  - 5µl biotin mix\* or 0.5ul dUTP-Cy3 1mM NOTE
  - 10µl 0.1M beta-mercaptoethanol
  - dilute (immediately before use) 1µl of DNaseI (2Uul)in 1ml distilled water;  
add 20µl to the nick-translation mix:  
the DNaseI should be calibrated to give fragments of 100-500bp
  - 10U DNA polymerase I
  - distilled sterile water to 100µl

2. incubate at 15°C for 2h
3. place at 4°C until it is checked on gel.

4. Take 5µl aliquot of each sample, add 4µl H<sub>2</sub>O and 1µl of 10x loading buffer. Run on 1% agarose gel. Load also a suitable MW marker (as PhiX-HaeIII). Inspect the gel on UV transilluminator: fragments should be between 100 and 500bp. If fragments are larger, add additional DNaseI to the samples, incubate for additional 30' and check again.
5. Stop reaction adding 4µl 0.5M EDTA

\* nick translation buffer:

10x buffer:

0.5M Tris-HCl pH 7.8-8

50mM MgCl<sub>2</sub>

0.5mg/ml BSA

\* dNTPs mix:

0.5mM dATP

0.5mM dCTP

0.5mM dGTP

\* biotin-16-dUTP mix

0.5mM dTTP

0.5mM bio-16-dUTP (Boehringer Mannheim)

NOTE: Direct labeling with dUTP-Cy3 (Amersham) can be used as an efficient alternative to biotin labeling. If co-hybridization are performed we usually label the first probe with dUTP-Cy3 and the second probe with biotin (detected with avidin-FITC).

---

## Fluorescence in situ hybridization (FISH)

### Slide denaturation

1. Prepare 200ul/slide of denaturing solution (70% deionized formamide/2xSSC).
2. Pre-warm slides at 60°C in dry oven.
3. Put 200ul of denaturation solution on each slide and incubate for exactly 2 min at 80°C in a dry oven or on a thermoblock plate.
4. Dehydrate slides in 70%, 90% and 100% ethanol, 3min each time (70% ethanol at -20°C).
5. Dry slides after dehydration.

### Probe denaturation

1. Precipitate 20µl labelled DNA (400ng) with 5µg human COT1 DNA (BRL), 3µg salmon sperm DNA, 2.8µl 3M Na acetate and 3 vol cold (-20°C) ethanol. Leave at -80°C for 15 minutes. Spin 15' (14,000RPM) at 4°C. Dry **completely** the pellet on a

- Savant fuge for few minutes.
2. Prepare hybridization mix (10  $\mu$ l per slide). Add to a test tube: 5 $\mu$ l of deionized formamide, 2 $\mu$ l dextran sulphate (50% in distilled water, autoclaved), 2 $\mu$ l distilled water and 1 $\mu$ l 20xSSC. If more slides have to be hybridized, a master mix can be prepared.
  3. Resuspend pellet in 10 $\mu$ l hybridization mix, by vortexing accurately;
  4. Denature DNA mix at 80°C for 8 min; transfer to 37°C for 20 min. Place on ice until used.
  5. NOTE: the amount of probe depends on the type of probe: 20-50 ug for repetitive DNA; 500-600 for painting probes.

## Hybridization

1. Apply 10 $\mu$ l hybridization mix to denatured slides, avoiding air bubbles.
2. Cover with 24x24mm clean coverslip; seal with rubber cement.
3. Incubate in a moist chamber at 37°C overnight.

## Post-hybridization washing and detection.

Do not allow slides to dry at any passages! All washings are done in Choplin jar.

1. Remove coverslips and wash 3 times for 5min in prewarmed solution (50% form. / 2xSSC) in a Choplin jar in a shaking waterbath at 42°C (we usually skip this first washing) ✓
2. Wash 3 times for 5' in prewarmed 0.1xSSC in waterbath at 60°C, **then go to 6 if direct labeling has been performed.** ✓
3. Apply 200 $\mu$ l of blocking solution per slide (3% BSA/4xSSC/0.1 Tween 20 ); cover with 24x60mm coverslip; transfer the slides in a moist chamber; incubate for 30' at 37°C.
4. Dilute stock solution of avidin-Cy3 in detection buffer (1%BSA/1xSSC/0.1 Tween 20) (Cy3 is from Amersham; Cy3 is stronger and more stable than other fluorochromes; use the same filter for rhodamine, or a specific one). Let coverslips slide off; apply 200 $\mu$ l detection solution per slide. Cover with 24x60mm coverslips. Transfer the slides in a dark moist chamber. Incubate at 37°C for 30 min.
5. Remove the coverslips; rinse the slides 3 times for 5 min in pre-warmed washing solution (4xSSC / 0.1 Tween 20) in waterbath at 42°C.
6. Counterstain with DAPI (200ng/ml in 2xSSC) or with propidium iodide (200ng/ml in 2xSSC). ✓
7. Rinse 2 min in 2xSSC 0.05% Tween 20 at room temperature.
8. Apply 30 $\mu$ l of antifade-mounting medium\* and cover with 24x60mm coverslip; slides can be stored for weeks in the dark at 4°C. ✓

\* Antifade-mounting medium (for 10ml):

- 0.233gr of DABCO (1,4-diazabicyclo-(2.2.2)octane, Sigma)
- 800 $\mu$ l H<sub>2</sub>O
- 200 $\mu$ l 1M Tris-HCl
- 9ml glycerol

Anandi Ravinath  
Asst. Manager  
Biochemicals Dept.



*Biochemicals*

**BOEHRINGER MANNHEIM** India Limited

54-A, Mathuradas Vasanji Road  
Chakala, Andheri (East)  
Bombay 400 093  
India

Tel: +91 (22) 837 07 94  
Telex: 011-79141 BM IN  
Fax: +91 (22) 837 99 06

Kashnaja, A.P.  
1.8.2000

**FISH protocol** 1st day  
for two-three chromosomes with pan-centromeric probe

1

**Pretreatment:**

- |      |   |      |           |   |
|------|---|------|-----------|---|
| ✓ 1  | Wash slides   | PBS  | 5 min     | RT  |
| ✓ 2  | Alcohol   | 70%  | 2-5 min   | RT  |
|      |   | 90%  | 2-5 min   | RT  |
|      |   | 100% | 2-5 min   | RT  |
| 3    | Air dry   |      |           | RT  |
| 4    | Pipette 100 ul of RNase A (100 ug/ml/2xSSC) per slide, overlay with coverslip   |      |           |   |
|      | Incubate in the moist chamber   |      | 60 min    | 37oC  |
|      | For 5 slides:   |      |           |   |
|      | { 445 ul H <sub>2</sub> O   |      |           |   |
|      | { 50 ul 20xSSC  |      |           |   |
|      | { 5 ul RNase 10 mg/ml   |      |           |   |
|      |   |      |           | <i>prepare in advance<br/>keep at -20°C</i> |
| ✓ 5  | 2xSSC   |      | 3 x 5 min | RT  |
|      | During 1st washing remove coverslips  |      |           |   |
|      | For 3 washes I usually prepared 400ml (40 ml 20xSSC + 360 ml water)   |      |           |   |
| ✓ 6  | PBS   |      | 5 min     | RT  |
| 7    | Pepsin (0.005% / 10mM HCl)  |      | 8-10 min  | 37oC  |
|      | { 99 ml H <sub>2</sub> O  |      |           |   |
|      | { 1 ml 1N HCl   |      |           |   |
|      | { 50 ul 10% pepsin  |      |           |   |
|      |   |      |           | <i>Prepare in advance<br/>keep at -20°C</i> |
|      |   |      |           | <i>44.5 ml<br/>0.5 ml<br/>2.5 ml</i>        |
|      | <b>NB:</b> solution should be pre-warmed in water-bath at least 30 min before use;<br>I usually prepared 100 ml 10mM HCL, pre-warmed it and added pepsin 5-10 min before use;<br>(if pepsin will stay for longer time at 37oC - it can lose activity) |      |           |   |
| 8    | PBS   |      | 5 min     | RT  |
| 9    | 50mM MgCl <sub>2</sub> / PBS  |      | 5 min     | RT  |
|      | { 5 ml 1M MgCl <sub>2</sub>   |      |           |   |
|      | { 95 ml PBS   |      |           |   |
| 10   | 1% formaldehyde/PBS/MgCl <sub>2</sub>   |      | 10 min    | RT  |
|      | { 3 ml 37% formaldehyde (commercial)  |      |           |   |
|      | { 5 ml 1M MgCl <sub>2</sub>   |      |           |   |
|      | { 92 ml PBS   |      |           |   |
|      |   |      |           | <i>twice 5 min RT</i>                       |
| ✓ 11 | PBS   |      | 2 x 5 min | RT  |
| ✓ 12 | Alcohol   | 70%  | 2-5 min   | RT  |
|      |   | 90%  | 2-5 min   | RT  |
|      |   | 100% | 2-5 min   | RT  |
|      | Air dry   |      |           | RT  |

*50 ul RNase.  
222.5 ul H<sub>2</sub>O  
25 ul 20xSSC  
2.5 ul RNase.  
10 mg/ml*

**NB:** for steps - 1, 2, 5, 6, 11, 12 - use shaker

**NB:** steps - 1, 2, 5, 6, 7, 8, 9, 10, 11, 12 - slides are washing in solutions in coplin jar



**Prehybridization of slides:**

(start 30 min before the end of paint's competition)

- 1 Prepare in Eppendorf tube (for 5 slides)

{ 350 ul deionize 100% formamide  
50 ul 0.5 M phosphate buffer (PB) *50mM*  
50 ul 20xSSC  
50 ul H<sub>2</sub>O

- 2 Pipette 100 ul of 70% formamide/2xSSC/50mM phosphate buffer per slide and overlay with coverslip 24 x 60 mm
- 3 Denature slides at 80°C for 2.5-3 min on the hot plate
- 4 Wash slides with ethanol series (during 1st wash remove coverslip)

70% ethanol	5 min (should be kept at -20°C before use) ✓
90% ethanol	2-5 min RT
100% ethanol	2-5 min RT
- 5 air dry

**Hybridization**

- 4 Mix all paints and CP in one Eppendorf tube mix very well, spin down for few seconds (several times !!!)
- 5 Pipette 20 ul of the mix solution onto each slide, overlay with coverslip 24 x 50 mm seal with Cow Gum, give time to glue to dry
- 6 Incubate overnight at 42°C in the moist chamber (incubation could be 2 overnights, if necessary)

FISH protocol

Detection:

I Make solutions for washing:

1 50% formamide / 2xSSC for 3 washing

200 ml 100% formamide + 20 ml 20xSSC + 180 ml water  
adjust pH 7.0

*Hand  
160ml*

2 0.1xSSC for 3 washing

3 0.05% Tween / 4xSSC (2-2,5 L)

{	200 ml	20xSSC
	5 ml	10% Tween 20
	800 ml	H2O
	<hr/> 1 L	

*0.1% Tween / 2xSSC*

<i>100ml</i>	<i>- 20xSSC</i>	<i>- 10ml</i>
<i>10ml</i>	<i>-</i>	<i>- 1ml</i>
<i>890ml</i>	<i>- H2O</i>	<i>- 89ml</i>

Pre-warm all solutions to 42-45°C in water bath

II Make antibodies for amplification:

Dilute Blocking Protein (BP from Cambio) to 10-15% (v/v) in wash solution from step 3.

This is solution A.

NB: sometimes 10-15% is too strong and could give big background, you can go down to 7%

Use solution A for dilution of antibodies:

I layer B 3 1 : 500 1 ul B 3 in 500 ul A

*( B 3 - Texas-Red Avidin )*

II layer B 4 1 : 250 { 2 ul B 4 and 2.5 ul F 1  
F 1 1 : 200 in 500 ul A

B 4 - Biotinylated goat anti-Avidin  
F 1 - Rabbit anti-FITC

III layer B 3 1 : 500 { 1 ul B 3 and 5 ul F 2  
F 2 1 : 100 in 500 ul A

F 2 - FITC goat anti-rabbit IgG

Incubate in the dark for 10 min at RT

Microfuge at 11,000g for 10 min and use supernatants

*amplifak*

*( 150ml ) 15ml - 1.2% Tween  
3ml -  
132 150ml )*

**Stringency washes:**

Carefully removed dried Cow Gum and wash slides

1 50% formamide / 2xSSC      3 x 5 min      42oC (45oC - Optional)  
During 1st washing remove coverslips

2 0.1xSSC      3 x 5 min      42oC (60oC - Optional)

3 0.05% Tween / 4xSSC      5 min      42oC

4 Pipette 100 ul of solution A per slide and overlay with coverslip 24x50 mm  
Incubate in moist chamber for 15-20 min at 37oC

NB: incubation can be up to 30 min

5 Wash slides  
0.05% Tween / 4xSSC      2 x 5 min      42oC

6 Pipette 100 ul of I layer of antibody per slide and overlay with coverslip  
Incubate in moist chamber for 20-30 min at 37oC

7 Wash slides  
0.05% Tween / 4xSSC      3 x 5 min      42oC

8 Pipette 100 ul of II layer of antibodies per slide and overlay with coverslip  
Incubate in moist chamber for 20-30 min at 37oC

9 Wash slides  
0.05% Tween / 4xSSC      3 x 5 min      42oC

10 Pipette 100 ul of III layer of antibodies per slide and overlay with coverslip  
Incubate in moist chamber for 20-30 min at 37oC

11 Wash slides  
0.05% Tween / 4xSSC      3 x 5 min      42oC

Repeat steps 8-11

Alcohol 70%      2-5 min      RT

90%      2-5 min      RT

100%      2-5 min      RT

Air dry      RT

Counterstain with DAPI (0.15 ug/ml in Vectashield)  
25 ul per slide under coverslip

Cover with coverslip

**Note:** Do not allow slides to dry out at any stage until last step - Air dry

Steps 5, 7, 9, 11 - washing with buffer - during 1st washing remove coverslip

1  
Bleaching

2

counterstaining

## Solutions

6

Prepare all solutions with double-distilled, deionized water

**PBS** (phosphate-buffer saline)

for 10 L:	NaCl	80 g	
	Na <sub>2</sub> HPO <sub>4</sub> · 2H <sub>2</sub> O	12 g	8 g
	KCl	2g	1.2
	KH <sub>2</sub> PO <sub>4</sub>	2g	.2

add water till 10 L

mix on the magnetic stir (approximately 1/2 hour)

adjust pH 7.2-7.4 (with 1M HCl)

put in 1L bottles and sterilize

store near 1 month at room temperature (RT)

<b>20 x SSC</b>	NaCl	175.3 g
	sodium citrate	88.24 g

add 950 ml water

mix on the magnetic stir (approximately 1/2 hour)

adjust pH 7.0 (with 1M or 5 M HCl)

add water till 1 L and sterilize

store near 1 month at RT

<b>1 M MgCl<sub>2</sub></b>	M.v. 203.3	
	for 250 ml MgCl <sub>2</sub>	<u>50.82g</u>

add water till 250 ml

mix on the magnetic stir

sterilize

store near 1 month at RT

<b>10% Tween</b>	100% Tween 20	10 ml
	water	90 ml

mix well

sterilize by filtration through a 0.2-micro m filter

store near 1 month at RT

Fluorescent antibody enhancer for DIG detection  
 10x PBS. pH 7.4.  
 1L. 80g NaCl  
 2g KCl.  
 12.4g Na<sub>2</sub>HPO<sub>4</sub>.  
 2.4g KH<sub>2</sub>PO<sub>4</sub>. (Boehringer)

**0.5 M PB** (phosphate buffer)prepare solution 1  
solution 20.5 M  $\text{Na}_2\text{HPO}_4$ 0.5 M  $\text{NaH}_2\text{PO}_4$ 

solution 1

 $\text{Na}_2\text{HPO}_4$ 

7.1 g in 100 ml water

should be dissolved on the magnetic stir  
pH will be approximately 9

solution 2

 $\text{NaH}_2\text{PO}_4 \cdot 1 \text{ H}_2\text{O}$   
pH will be approximately 4

6.9 g in 100 ml water

in full volume of solution 1 (which is still on the magnetic stir)  
add solution 2 to adjust pH till 7.0  
(when I did it, for 100 ml solution 1 approximately - 40 ml of solution 2)  
dispense into 1 ml aliquots and store at -20°C

**100% deionized formamide**

mix 50 ml of 100% (commercial) formamide and 5 g Serdolit MB  
stir on the magnetic stir - 1,5-2 hours at RT  
filter twice through Whatman No.1 filter paper  
dispense into 1 ml aliquots and store at -20°C

Serdolit MB - ion-exchange resin; SERVA (company), cat.number 45500

**10% pepsin**

1 g pepsin dissolve in 10 ml of 10 mM HCl  
mix well on vortex  
dispense into 1 ml aliquots and store at -20°C

Krishna, A.P.  
1-8-2000  
(J)

**FISH protocol** 1st day  
for two-three chromosomes with pan-centromeric probe

**Pretreatment:**

1	Wash slides	PBS	5 min	RT
2	Alcohol	70%	2-5 min	RT
		90%	2-5 min	RT
		100%	2-5 min	RT
3	Air dry			RT
4	Pipette 100 ul of RNase A (100 ug/ml/2xSSC) per slide, overlay with coverslip			
	Incubate in the moist chamber		60 min	37oC
	For 5 slides:			
	445 ul H <sub>2</sub> O			
	50 ul 20xSSC			
	5 ul RNase 10 mg/ml			
5	2xSSC		3 x 5 min	RT
	During 1st washing remove coverslips			
	For 3 washes I usually prepared 400ml (40 ml 20xSSC + 360 ml water)			
6	PBS		5 min	RT
7	Pepsin (0.005% / 10mM HCl)		8-10 min	37oC
	99 ml H <sub>2</sub> O			
	1 ml 1N HCl			
	50 ul 10% pepsin			
	<b>NB:</b> solution should be pre-warmed in water-bath at least 30 min before use;			
	I usually prepared 100 ml 10mM HCL, pre-warmed it and added pepsin 5-10 min before use;			
	(if pepsin will stay for longer time at 37oC - it can lose activity)			
8	PBS		5 min	RT
9	50mM MgCl <sub>2</sub> / PBS		5 min	RT
	5 ml 1M MgCl <sub>2</sub>			
	95 ml PBS			
10	1% formaldehyde/PBS/MgCl <sub>2</sub>		10 min	RT
	3 ml 37% formaldehyde (commercial)			
	5 ml 1M MgCl <sub>2</sub>			
	92 ml PBS			
11	PBS		2 x 5 min	RT
12	Alcohol	70%	2-5 min	RT
		90%	2-5 min	RT
		100%	2-5 min	RT
	Air dry			RT

**NB:** for steps - 1, 2, 5, 6, 11, 12 - use shaker

**NB:** steps - 1, 2, 5, 6, 7, 8, 9, 10, 11, 12 - slides are washing in solutions in coplin jar



**Prehybridization of slides:**

(start 30 min before the end of paint's competition)

- 1 Prepare in Eppendorf tube (for 5 slides)  
  
350 ul deionize 100% formamide  
50 ul 0.5 M phosphate buffer (PB)  
50 ul 20xSSC  
50 ul H<sub>2</sub>O
- 2 Pipette 100 ul of 70% formamide/2xSSC/50mM phosphate buffer per slide and overlay with coverslip 24 x 60 mm
- 3 Denature slides at 80oC for 2.5-3 min on the hot plate
- 4 Wash slides with ethanol series (during 1st wash remove coverslip)  
  
70% ethanol                      5 min (should be kept at -20oC before use)  
90% ethanol                      2-5 min    RT  
100% ethanol                      2-5 min    RT
- 5 air dry

**Hybridization**

- 4 Mix all paints and CP in one Eppendorf tube  
mix very well, spin down for few seconds (several times !!!)
- 5 Pipette 20 ul of the mix solution onto each slide,  
overlay with coverslip 24 x 50 mm  
seal with Cow Gum, give time to glue to dry
- 6 Incubate overnight at 42oC in the moist chamber  
(incubation could be 2 overnights, if necessary)

## FISH protocol

4

### Detection:

#### I *Make solutions for washing:*

- 1 50% formamide / 2xSSC for 3 washing  
200 ml 100% formamide + 20 ml 20xSSC + 180 ml water  
adjust pH 7.0
- 2 0.1xSSC for 3 washing
- 3 0.05% Tween / 4xSSC (2-2,5 L)  

200 ml	20xSSC
5 ml	10% Tween 20
800 ml	H2O
<hr style="width: 100%;"/>	
1 L	

Pre-warm all solutions to 42-45°C in water bath

#### II *Make antibodies for amplification:*

Dilute Blocking Protein (BP from Cambio) to 10-15% (v/v) in wash solution from step 3.  
This is solution A.

NB: sometimes 10-15% is too strong and could give big background, you can go down to 7%

Use solution A for dilution of antibodies:

I layer	B 3	1 : 500	1 ul B 3 in 500 ul A	B 3 - Texas-Red Avidin
II layer	B 4	1 : 250	2 ul B 4 and 2.5 ul F 1 in 500 ul A	B 4 - Biotinylated goat anti-Avidin
	F 1	1 : 200		F 1 - Rabbit anti-FITC
III layer	B 3	1 : 500	1 ul B 3 and 5 ul F 2 in 500 ul A	F 2 - FITC goat anti-rabbit IgG
	F 2	1 : 100		

Incubate in the dark for 10 min at RT

Microfuge at 11,000g for 10 min and use supernatants

**Stringency washes:**

Carefully removed dried Cow Gum and wash slides

- 1 50% formamide / 2xSSC            3 x 5 min            42oC (45oC - Optional)  
During 1st washing remove coverslips
- 2 0.1xSSC                                3 x 5 min            42oC (60oC - Optional)
- 3 0.05% Tween / 4xSSC            5 min                42oC
- 4 Pipette 100 ul of solution A per slide and overlay with coverslip 24x50 mm  
Incubate in moist chamber for 15-20 min at 37oC
- NB: incubation can be up to 30 min
- 5 Wash slides  
0.05% Tween / 4xSSC            2 x 5 min            42oC
- 6 Pipette 100 ul of I layer of antibody per slide and overlay with coverslip  
Incubate in moist chamber for 20-30 min at 37oC
- 7 Wash slides  
0.05% Tween / 4xSSC            3 x 5 min            42oC
- 8 Pipette 100 ul of II layer of antibodies per slide and overlay with coverslip  
Incubate in moist chamber for 20-30 min at 37oC
- 9 Wash slides  
0.05% Tween / 4xSSC            3 x 5 min            42oC
- 10 Pipette 100 ul of III layer of antibodies per slide and overlay with coverslip  
Incubate in moist chamber for 20-30 min at 37oC
- 11 Wash slides  
0.05% Tween / 4xSSC            3 x 5 min            42oC

Repeat steps 8-11

Alcohol	70%	2-5 min	RT
	90%	2-5 min	RT
	100%	2-5 min	RT
Air dry			RT

Counterstain with DAPI (0.15 ug/ml in Vectashield)  
25 ul per slide under coverslip

Cover with coverslip

**Note:** Do not allow slides to dry out at any stage until last step - Air dry

Steps 5, 7, 9, 11 - washing with buffer - during 1st washing remove coverslip

Prepare all solutions with double-distilled, deionized water

**PBS** (phosphate-buffer saline)

for 10 L:	NaCl	80 g
	Na <sub>2</sub> HPO <sub>4</sub> · 2H <sub>2</sub> O	12 g
	KCl	2g
	KH <sub>2</sub> PO <sub>4</sub>	2g
	<hr/>	
	add water till 10 L	

mix on the magnetic stir (approximately 1/2 hour)  
 adjust pH 7.2-7.4 (with 1M HCl)  
 put in 1L bottles and sterilize  
 store near 1 month at room temperature (RT)

<b>20 x SSC</b>	NaCl	175.3 g
	sodium citrate	88.24 g
	<hr/>	
	add 950 ml water	

mix on the magnetic stir (approximately 1/2 hour)  
 adjust pH 7.0 (with 1M or 5 M HCl)  
 add water till 1 L and sterilize  
 store near 1 month at RT

<b>1 M MgCl<sub>2</sub></b>	M.v. 203.3	
	for 250 ml MgCl <sub>2</sub>	50.82g
	add water till 250 ml	
	mix on the magnetic stir	
	sterilize	
	store near 1 month at RT	

<b>10% Tween</b>	100% Tween 20	10 ml
	water	90 ml
	mix well	
	sterilize by filtration through a 0.2-micro m filter	
	store near 1 month at RT	

**0.5 M PB (phosphate buffer)**

prepare      solution 1                      0.5 M  $\text{Na}_2\text{HPO}_4$   
                  solution 2                      0.5 M  $\text{NaH}_2\text{PO}_4$

solution 1                                       $\text{Na}_2\text{HPO}_4$                       7.1 g in 100 ml water  
should be dissolved on the magnetic stir  
pH will be approximately 9

solution 2                                       $\text{NaH}_2\text{PO}_4 \cdot 1 \text{ H}_2\text{O}$                       6.9 g in 100 ml water  
pH will be approximately 4

in full volume of solution 1 (which is still on the magnetic stir)  
add solution 2 to adjust pH till 7.0  
(when I did it, for 100 ml solution 1 approximately - 40 ml of solution 2)  
dispense into 1 ml aliquots and store at  $-20^\circ\text{C}$

**100% deionized formamide**

mix 50 ml of 100% (commercial) formamide and 5 g Serdolit MB  
stir on the magnetic stir - 1,5-2 hours at RT  
filter twice through Whatman No.1 filter paper  
dispense into 1 ml aliquots and store at  $-20^\circ\text{C}$

Serdolit MB - ion-exchange resin; SERVA (company), cat.number 45500

**10% pepsin**

1 g pepsin dissolve in 10 ml of 10 mM HCl  
mix well on vortex  
dispense into 1 ml aliquots and store at  $-20^\circ\text{C}$

Krishna, A.P.  
1-8-2000  
(J)

**FISH protocol**                      **1st day**  
short version for one chromosome and pan-centromeric probe

**Pretreatment:**

* PBS	5 min	RT
* 70% ethanol	5 min	RT
90% ethanol	5 min	RT
100% ethanol	5 min	RT
* Air dry		RT
* 100 ul of RNase A (100 ug /ml / 2xSSC) (stock RNase 10 mg / ml)	60 min	37oC
* 2xSSC	3 x 5 min	RT
* PBS	5 min	RT
* Pepsin (0.005% / 10mM HCl) (stock 10% pepsin)	10 min	37oC
* PBS	5 min	RT
* 50mM MgCl <sub>2</sub> / PBS	5 min	RT
* 1% formaldehyde/PBS/MgCl <sub>2</sub>	10 min	RT
* PBS	5 min	RT
* 70% ethanol	5 min	RT
90% ethanol	5 min	RT
100% ethanol	5 min	RT
* Air dry		RT

**Probe denaturation and competition:****A** Paint for specific chromosome

warm paint to 42°C

mix 3 ul of paint with 7 ul of hybridization buffer (HB)

denature paint at 65°C for 10 min

incubate for 60 to 90 min at 37°C

**B** Pan-centromeric probe (CP)

warm CP and HB at 37°C for 5-10 min

mix 2 ul of CP with 7 ul of HB

denature CP at 85°C for 10 min

put on ice for 2-3 min - then spin down for few seconds and keep on ice

**Prehybridization of slides:**

- |   |  |       |                            |
|---|--|-------|----------------------------|
| 1 | 100 ul/slide of denaturation solution<br>(70% formamide / 2xSSC / 50mM phosphate buffer) | 80°C  | 2.5-3 min                  |
| 2 | 70% ethanol  | 5 min | (kept at -20°C before use) |
|   | 90% ethanol  | 5 min | RT                         |
|   | 100% ethanol   | 5 min | RT                         |
| 3 | Air dry  |       |                            |

**Hybridization**

- 4 Mix paint and CP in one Eppendorf tube
- 5 Pipette 20 ul of the mix solution on slide, overlay with coverslip, seal with Cow Gum
- 6 Incubate overnight at 42°C in the moist chamber

**FISH protocol****2nd day****3****Detection:**

- 1 Make a wash solution 0.05% Tween / 4xSSC
- 2 Dilute Blocking Protein (BP) to 10-15% (v/v) in wash solution  
This is solution A.
- 3 Use solution A for dilution of antibodies:
 

I layer	B 3	1 : 500	B 3 - Texas-Red Avidin
II layer	B 4	1 : 250	B 4 - Biotinylated goat anti-Avidin
	F 1	1 : 200	F 1 - Rabbit anti-FITC
III layer	B 3	1 : 500	F 2 - FITC goat anti-rabbit IgG
	F 2	1 : 100	

Keep antibodies in the dark for 10 min at RT

Microfuge at 11,000g for 10 min and use supernatants

- 4 Pre-warm to 45oC:

50% formamide / 1xSSC, pH 7.0  
0.1xSSCfor 3 washing  
for 3 washing

- |    |   |           |      |
|----|---|-----------|------|
| 5  | 50% formamide / 1xSSC                     | 2 x 5 min | 45oC |
| 6  | 0.1xSSC                                   | 3 x 5 min | 45oC |
| 7  | 0.05% Tween / 4xSSC                       | 5 min     | 45oC |
| 8  | 100 ul / slide of solution A              | 20 min    | 37oC |
| 9  | 0.05% Tween / 4xSSC                       | 2 x 5 min | 45oC |
| 10 | 100 ul / slide of I layer of antibody     | 20 min    | 37oC |
| 11 | 0.05% Tween / 4xSSC                       | 3 x 5 min | 45oC |
| 12 | 100 ul / slide of II layer of antibodies  | 20 min    | 37oC |
| 13 | 0.05% Tween / 4xSSC                       | 3 x 5 min | 45oC |
| 14 | 100 ul / slide of III layer of antibodies | 20 min    | 37oC |
| 15 | 0.05% Tween / 4xSSC                       | 3 x 5 min | 45oC |

Repeat steps 12-15

- |    |                        |      |    |
|----|------------------------|------|----|
| 16 | Alcohol 70%            | 5min | RT |
|    | 90%                    | 5min | RT |
|    | 100%                   | 5min | RT |
| 17 | Air dry                |      | RT |
| 18 | Counterstain with DAPI |      |    |

Krishna's

25-4-2000

<sup>1</sup>  
(From  
Janna)

## FISH protocol

### Pretreatment:

1	Wash slides	PBS	5 min	RT
2	Alcohol	70%	2-5 min	RT
		90%	2-5 min	RT
		100%	2-5 min	RT
3	Air dry			RT
4	Pipette 100 $\mu$ l of RNase A (100 $\mu$ g/ml) per slide, overlay with coverslip			
	Incubate in the moist chamber		60 min	37oC
	{ 445 $\mu$ l H <sub>2</sub> O			
	50 $\mu$ l 20xSSC			
	5 $\mu$ l RNase 10 mg/ml (store at -20oC)			
5	2xSSC		3 x 5 min	RT
	During 1st washing remove coverslips			
6	PBS		5 min	RT
7	Pepsin (0.005% / 10mM HCl)		8-10 min	37oC
	{ 99 ml H <sub>2</sub> O			
	1 ml 1N HCl			
	50 $\mu$ l 10% pepsin (store at -20oC)			
<b>NB:</b> solution should be pre-warmed in water-bath before use at least 30 min before use; I usually prepare 10mM HCl, pre-warmed it and add pepsin 10 min before putting slides in solution (if pepsin will stay for long time at 37oC - it can lose activity)				
8	PBS		5 min	RT
9	50mM MgCl <sub>2</sub> / PBS		5 min	RT
	{ 5 ml 1M MgCl <sub>2</sub>			
	95 ml PBS			
10	1% formaldehyde/PBS/MgCl <sub>2</sub>		10 min	RT
	{ 3 ml 37% formaldehyde			
	5 ml 1M MgCl <sub>2</sub>			
	92 ml PBS			
11	PBS		2 x 5 min	RT
12	Alcohol	70%	2-5 min	RT
		90%	2-5 min	RT
		100%	2-5 min	RT
	Air dry			RT

**NB:** for steps - 1, 2, 5, 6, 11, 12 - use shaker



**Prehybridization of slides:**

(start 30 min before the end of paint's competition)

- 1 Pipette 100  $\mu$ l of 70% formamide/2xSSC/50mM phosphate buffer per slide and overlay with coverslip
  - { 350  $\mu$ l deionize 100% formamide (store at -20oC)
  - 50  $\mu$ l 0.5 M phosphate buffer (store at -20oC)
  - 50  $\mu$ l 20xSSC
  - 50  $\mu$ l H<sub>2</sub>O
- 2 Denature slides at 80oC for 2.5-3 min on the hot plate
- 3 

70% ethanol	5 min	should be kept at -20oC before use
90% ethanol	2-5 min	RT
100% ethanol	2-5 min	RT

air dry

**Hybridization**

- 4 Mix all paints and CP in one Eppendorf tube  
**Note:** mix very well (several times), spin down for few seconds
- 5 Pipette 20  $\mu$ l of the mix solution onto each slide, overlay with coverslips, seal with Cow Gum, give time to glue to dry
- 6 Incubate overnight at 42oC in the moist chamber (incubation could be 2 overnights, if necessary)

**Detection:**

## 1 Make a wash solution:

a. add 0.5 ml of detergent to 1L of 4xSSC (Cambio)

b. 0.05% Tween / 4xSSC

{	200 ml 20xSSC
	5 ml 10% Tween 20
	800 ml H <sub>2</sub> O

NB: I am usually use - b.

## 2 Dilute Blocking Protein (BP) to 10-15% (v/v) in wash solution from step 1.

This is solution A.

NB: sometimes 10-15% is too strong and give big background, you can go down to 7%

## 3 Use solution A for dilution of antibodies:

I layer B 3 1 : 500

B 3 - Texas-Red Avidin

II layer B 4 1 : 250

B 4 - Biotinylated goat anti-Avidin

F 1 1 : 200

F 1 - Rabbit anti-FITC

III layer B 3 1 : 500

F 2 - FITC goat anti-rabbit IgG

F 2 1 : 100

NB: 100 µl per hybridisation area plus 50 µl excess

## 4 Incubate in the dark for 10 min at RT

Microfuge at 11,000g for 10 min and use supernatants

## 5 Pre-warm to 42-45°C:

2xSSC

for 1 washing

50% formamide / 2xSSC, pH 7.0

for 2-3 washing

0.1xSSC

for 3 washing

## 6 Carefully removed dried Cow Gum and rinse off coverslips in the jar of 2xSSC

**Note:** Do not allow slides to dry out at any stage until step 16

**Stringency washes:**

- 7 Wash slides
- |                       |             |      |                   |
|-----------------------|-------------|------|-------------------|
| 50% formamide / 2xSSC | 2-3 x 5 min | 42oC | (45oC - Optional) |
| 0.1xSSC               | 3 x 5 min   | 42oC | (60oC - Optional) |
| 0.05% Tween / 4xSSC   | 5 min       | 42oC |                   |
- 8 Pipette 100  $\mu$ l of solution A per slide and overlay with coverslip  
Incubate in moist chamber for 15-20 min at 37oC
- NB: incubation can be up to 30 min
- 9 Wash slides
- |                     |           |      |  |
|---------------------|-----------|------|--|
| 0.05% Tween / 4xSSC | 2 x 5 min | 42oC |  |
|---------------------|-----------|------|--|
- 10 Pipette 100  $\mu$ l of I layer of antibody per slide and overlay with coverslip  
Incubate in moist chamber for 20-30 min at 37oC
- 11 Wash slides
- |                     |           |      |  |
|---------------------|-----------|------|--|
| 0.05% Tween / 4xSSC | 3 x 5 min | 42oC |  |
|---------------------|-----------|------|--|
- 12 Pipette 100  $\mu$ l of II layer of antibodies per slide and overlay with coverslip  
Incubate in moist chamber for 20-30 min at 37oC
- 13 Wash slides
- |                     |           |      |  |
|---------------------|-----------|------|--|
| 0.05% Tween / 4xSSC | 3 x 5 min | 42oC |  |
|---------------------|-----------|------|--|
- 14 Pipette 100  $\mu$ l of III layer of antibodies per slide and overlay with coverslip  
Incubate in moist chamber for 20-30 min at 37oC
- 15 Wash slides
- |                     |           |      |  |
|---------------------|-----------|------|--|
| 0.05% Tween / 4xSSC | 3 x 5 min | 42oC |  |
|---------------------|-----------|------|--|
- Repeat steps 12-15
- 16 Alcohol
- |      |         |    |
|------|---------|----|
| 70%  | 2-5 min | RT |
| 90%  | 2-5 min | RT |
| 100% | 2-5 min | RT |
- 17 Air dry RT
- 18 Counterstain with DAPI (0.15  $\mu$ g/ml in Vectashield)  
25  $\mu$ l per slide under coverslip
- 19 Cover with coverslip,

krishna's 25-4-2000

## Protocol A

**Approx. time: Slide Preparation: 30+overnight  
(or +90 min) + 30 min  
Probe preparation: 75 min**

# cambio

34 Newnham Road, Cambridge CB3 9EY  
Tel: +44 (0)1223 366500 Fax: +44 (0)1223 350069  
e-mail: support@cambio.demon.co.uk  
http://www.cambio.co.uk

## IN SITU HYBRIDISATION PROTOCOL

This protocol is for use with all human STAR★FISH chromosome paints.

Also for use with the Biotin (FITC) Painting Kit (1089-KB), Biotin (Texas Red) Painting Kit (1082-KT) and the FITC Amplification Kit (1084-KF).

### SLIDE PREPARATION

1. One day before needed, prepare fresh slides from 3:1 methanol:acetic acid fixed preparations. Air dry overnight. (Alternatively, prepare fresh slides from 3:1 methanol:acetic acid and age slides for 1.5 hours on a hot plate at 65°C).
2. Dehydrate by serial ethanol washing (70%(v/v), 70%, 90%, 90%, 100%) and air dry.
3. Denature slides by incubating in 70% (v/v) formamide (35ml formamide plus 15ml of 2X SSC) at 65°C for 1.5-2 min. Ensure that the 70% (v/v) formamide has reached temperature before use and do not exceed this time. For best results use fresh deionised formamide.
4. Quench in ice-cold 70% (v/v) ethanol and dehydrate using ethanol as in 2 above. Air dry.

### MEANWHILE

5. Warm paints to 42°C and mix well before use. Take sufficient probe (3µl of probe, plus hybridisation buffer to a total of 15µl per slide) and place in a small Eppendorf tube. Denature paint by incubating at 65°C for 10 mins then at 37°C for 60 to 90 mins. Note: Paints may be prediluted (1 part paint:4 parts hybridisation buffer) and stored in hybridisation buffer if so desired.
6. Pipette 15µl of probe onto each slide, overlay with cover slips, seal with Cow Gum and incubate overnight at 42°C.

### 1 X SSC Recipe:

0.015 M Na Citrate, pH 7.0 containing 0.15 M NaCl.

*Note: These products are not for use in humans and are for research purposes only.*



A QUALITY ASSURED ORGANIZATION

Revised 9th March 1998

Protocol F

*Approx. time: Slide Preparation: 30 min*

*Detection: 140 min*

## DETECTION PROTOCOL

### Dual Colour Painting Kit (Texas Red & FITC) (1090-KD)

for use with biotin and FITC labelled STAR★FISH chromosome paints.

1. Make a Wash Solution: add 250µl of Detergent to 500ml of 4X SSC.
2. Dilute Blocking Protein to 15% in the 4X SSC/Detergent Wash Solution from step 1. This is solution A.
3. Dilute detection reagents using Solution A:
  - a) Make a 1:500 dilution of Detection Reagent B3 (100µl per slide plus 50µl excess).
  - b) Make a mix of 1:250 final dilution of Detection Reagent B4 and a 1:200 final dilution of Detection Reagent F1 (100µl per slide plus 50µl excess).
  - c) Make a mix of 1:500 final dilution of Detection Reagent B3 and a 1:100 final dilution of Detection F2 (100µl per slide plus 50µl excess).
4. Incubate in the dark for 10 min at room temperature. Microfuge at 11,000g for 10 min and use supernatants.
5. Pre-warm to 45°C: one Coplin jar of 2X SSC, two Coplin jars of 50% (v/v) formamide/1X SSC (ie: 1:1 mix of formamide and 2X SSC) and two Coplin jars of 0.1X SSC.
6. Carefully remove dried Cow Gum and rinse off cover slips in the jar of 2X SSC for 5 mins.

NOTE: DO NOT ALLOW SLIDES TO DRY OUT AT ANY STAGE UNTIL STEP 16.

7. Stringency Washes: Wash slides by incubating for 5 min at 45°C in each of the two jars containing formamide.
8. Wash slides twice by incubating for 5 min at 45°C in each of the two jars containing 0.1X SSC.
9. Pipette 100µl of the diluted Detection Reagent B3 onto a slide and cover with a Parafilm cover slip. Incubate at 37°C for 15-20 min in a humidified box.
10. Wash 3 times in Wash Solution at 42°C for 5 min each time.
11. Pipette 100µl of diluted Detection Reagent B4/F1 supernatant onto each slide and overlay with Parafilm cover slip. Incubate in a humidified box for 15-20 min at 37°C.
12. Wash 3 times in Wash Solution at 42°C for 5 min each time.
13. Pipette 100µl of diluted Detection Reagent B3/F2 supernatant onto each slide and overlay with a Parafilm cover slip. Incubate in a humidified box for 15-20 min at 37°C.
14. Wash 3 times in wash solution at 42°C for 5 min each time.
15. Dilute 1µl of Counterstain 1 with 9µl of distilled water. Add 6µl of the diluted Counterstain 1 to 200µl of Mountant. Mix well. Discard unused diluted counterstain.
16. Briefly dehydrate slides using an ethanol series as in Step 3 (from slide preparation protocol), air dry and mount by adding 20µl of the Mountant/Counterstain mixture. Overlay with cover slip and seal with nail varnish.  
Store slides in the dark at 4°C.
17. View slides using standard epifluorescence filters for FITC and Texas Red and for DAPI.

1 X SSC Recipe:

0.015 M Na Citrate, pH 7.0 containing 0.15 M NaCl.

*Note: These products are not for use in humans and are for research purposes only.*



34 Newnham Road, Cambridge CB3 9EY  
Tel: +44 (0)1223 366500 Fax: +44 (0)1223 350069  
e-mail: [postmaster@cambio.demon.co.uk](mailto:postmaster@cambio.demon.co.uk)  
<http://www.cambio.co.uk>

**Protocol G**      **Approx. time : Probe Preparation: 20 min +  
16 hours (overnight) + 30 min**

### **IN SITU HYBRIDISATION PROTOCOL**

#### **STAR★FISH Pan Centromeric Chromosome Paint (1141).**

To be used with the Biotin (Texas Red) Painting Kit (1082-KT) or Dual Colour Painting Kit (1090-KD)

1. Warm probe to 37°C for 5 min and mix well.
2. Add 25µl of probe to a microcentrifuge tube.
3. Denature probe for 10 min at 85°C. Immediately chill on ice.
4. Apply 25µl of probe to slide.
5. Hybridise for 16 hours at 37°C.
6. **POST HYBRIDISATION WASH:**  
Remove cover slip in 2X SSC at 37°C. 5 min.  
Wash slides in 50% formamide/2X SSC 37°C, 2X 5 min.  
Wash slides in 2X SSC, 2X 5min.
7. Detect the probe using Cambio's Biotin (Texas Red) Painting Kit (1082-KT) or Dual Colour Painting Kit (1090-KD).

#### **1 X SSC RECIPE:**

0.015 M Na Citrate, pH 7.0 containing 0.05 M NaCl.

*Note: These products are not for use in humans and are for research purposes only.*



A QUALITY  
ASSURED ORGANIZATION